**DENDIA**  
DENTAL INSTRUMENTS

*...the first step  
to a perfect smile*

**DENTALINSTRUMENTE**

DENTAL INSTRUMENTS  
INSTRUMENTS DENTAIRE  
INSTRUMENTAL DENTAL

*Edition 2024*

*...the first step  
to a perfect smile*

## INTRODUCTION

EINLEITUNG | INTRODUCTION | INTRODUCCIÓN

002

**ABOUT DENDIA**

ÜBER DENDIA  
SOBRE DENDIA  
À PROPOS DE DENDIA

004

**GENERAL INFORMATION**

ALLGEMEINE INFORMATIONEN  
INFORMATIONS GÉNÉRALES  
INFORMACIÓN GENERAL

008

**CLEANING**

REINIGUNG  
NETTOYAGE  
LIMPIEZA

## PRAXIS

DENTISTRY | CABINET | CLINICA

010

**DICA DIAMOND**

DICA DIAMANT  
DIAMANTS DICA  
DIAMANTES DICA

064

**POLISHER**

POLIERER  
POLISSOIRS  
PULIDORES

042

**CARBIDE**

HARTMETALL  
CARBURE  
CARBURO

76

**ORTHO-DONTICS**

KIEFERORTHOPÄDIE  
ORTHODONTIE  
ORTODONTICA

052

**STEEL**

STAHL  
ACIER  
ACERO

080

**SURGERY**

CHIRURGIE  
CHIRURGIE  
CIRUGÍA

056

**ABRASIVES**

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KITS  
KITS  
KITS


170

**APPLICATION AND SAFETY INSTRUCTIONS**

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ÜBER DENDIA | SOBRE DENDIA | À PROPOS DE DENDIA

## DENDIA

**EN:** Since 1947 customers rely on the high-quality products made by Dendia. A high level of technical competence in connection with modern production methods and high-quality standards are the basis for the satisfaction of our customers.

Our wide product range containing all types of rotary instruments covers the needs of dentists and dental technicians all over the world.

Close and individually coordinated cooperation with our global partners as well as mutual trust are our keys to success on the market.

As a manufacturer of class IIa medical products Dendia of course is certified according to DIN EN ISO 13485:2016 and RL 93/42/EWG.

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
  
Dellagiacomma Ricardo

## EC Certificate

mdc medical device certification GmbH  
 Notified Body 0433  
 herewith certifies that

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for the scope  
 rotary instruments for dentistry applications  
 (see attachments)  
 has introduced and applies a  
**Quality System**  
 for the design, manufacture and final inspection.

The mdc audit has proven that this quality system  
 meets all requirements according to  
**Annex II – excluding Section 4**  
**of the Council Directive 93/42/EEC**  
 of 14 June 1993 concerning medical devices.  
 The surveillance will be held as specified in Annex II, Section 6

Valid from 2018-08-28  
 Valid until 2024-05-28  
 Registration no. D400900008  
 Report no. P18-02986-131237  
 Stuttgart 2018-08-28

*[Signature]*  
 Head of Certification Body


mdc medical device certification GmbH  
 Kriegerstraße 8  
 D-70191 Stuttgart, Germany  
 Phone: +49 (0)711 263697-0  
 Fax: +49(0)711-263697-10  
 Internet: <http://www.mdc-co.de>

mdc medical device certification GmbH  
 Kriegerstraße 8  
 D-70191 Stuttgart, Germany  
 Phone: +49(0)711-263697-0  
 Fax: +49(0)711-263697-10  
 Internet: <http://www.mdc-co.de>

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
→ *The recommended speed are a theoretical value, you can find the specific recommended speed in the product-related information.*

*To over exceed the given rotation speed can lead to impairment of the efficiency of the products and/or be of danger for the patient as well as the user. The choice of the working speed for a specific application case is determined by the material subsequently to be used, the instrument power, and the maximal down force pressure ultimately decided upon by the user.*

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







# CLEANING PROCEDURE OF MEDICAL PRODUCTS


REINIGUNG VON MEDIZINPRODUKTEN | NETTOYAGE DES PRODUITS MÉDICAUX |  
PROCEDIMIENTO DE LIMPIEZA DE LOS PRODUCTOS MÉDICOS

Newly delivered unsterile instruments must be prepared accordingly before first use!

1. Pre-treatment	
<ul style="list-style-type: none"> <li>› Abrasive impurities need to be removed from the products directly after use (within two hours maximum). To do so, use running water or a disinfectant solution (must not contain aldehydes, its effectiveness should be established, suitable and compatible for the products)</li> </ul>	
2a. Cleaning and disinfecting	
<ul style="list-style-type: none"> <li>› Disassemble the products</li> <li>› Place the disassembled products in the cleaning bath for the prescribed contact time and concentration, take care that the products do not touch each other</li> <li>› Products with lumen: rinse all instrument lumens 5 times at the beginning and or at the end of the contact time using a disposal syringe and a cannula</li> <li>› Take the products out of the cleaning bath and rinse at least 3 times thoroughly with water</li> <li>› Check the products</li> <li>› Dry the products by blowing them dry using filtered pressurised air (Use only filtered air for drying and water that is either sterile or low in germs and endotoxins)</li> <li>› Wrap the products</li> </ul>	
2b. Cleaning and disinfecting with disinfectant / RDG (MACHINE CLEANING)	
<ul style="list-style-type: none"> <li>› The effectiveness of the disinfectant has to be certified</li> <li>› Disassemble the products as far as possible</li> <li>› Place the disassembled products in the disinfectant and take care that the products do not touch each other</li> <li>› Start the program</li> <li>› Remove the products from the disinfectant after the program has finished</li> <li>› Check and wrap the products straight after removal if possible</li> <li>› Ensure that the disinfectant is regularly maintained and checked. Use only filtered air for drying and water that is either sterile or low in germs and endotoxins</li> </ul>	
3. Checking / Maintenance / Packaging	
<ul style="list-style-type: none"> <li>› Check the products for corrosion, damaged surfaces, bare patches, broken, chipped-off edges, deformations or impurities and eliminate damaged products. Products that are still contaminated need to be cleaned and disinfected once more</li> <li>› Re-assemble disassembled products, instrument oil must not be used</li> <li>› Arrange the cleaned and disinfected products in the dedicated bur block / sterilisation tray.</li> <li>› Wrap the products, bur blocks, sterilisation trays using disposable sterilisation packaging (disposable or double packaging) or sterilisation containers</li> </ul>	
4. Sterilisation	
<ul style="list-style-type: none"> <li>› Steam sterilisation: use of a fractional vacuum process or a gravitation process</li> <li>› Maximum sterilisation temperature 134°C (273°F)</li> <li>› At least 3 min (or 18 min at prion deactivation) at 132°C (270°F) / 134°C (273°F) and drying time of 20 min</li> </ul>	
5. Storage	
<ul style="list-style-type: none"> <li>› The products need to be stored in sterilisation wrapping in a dry and dust-free place. Please note the storage-life resulting from the validation of the sterilisation wrapping</li> </ul>	
<p><b>i</b></p> <ul style="list-style-type: none"> <li>› New unsteril supplied instruments must be prepared before first use.</li> <li>› Instruments made of tool steel are not suitable for sterilisation without an appropriate pretreatment.</li> <li>› When cleaning instruments unprotected against corrosion (like steel instruments) a corrosion-inhibiting disinfectants and cleaning agents must be used and they have to be pretreated with rust preventing spray before sterilisation.</li> <li>› Furthermore, make sure that instruments of different materials should never be reprocessed together.</li> <li>› Use particular care during the cleaning of the grinding surfaces of the instruments and ensure that all residues are removed by using a brush.</li> <li>› In case of severe contamination of the instrument, it is recommended to use an ultrasonic bath for cleaning.</li> <li>› Usage of protection gloves during work with contaminated instruments is highly recommended!</li> <li>› Bur blocks / instrument trays: <ul style="list-style-type: none"> <li>› Cleaning and disinfecting only without products being loaded (products must not be cleaned and disinfected whilst they are in the bur block / instrument tray)</li> <li>› Burblocks and other instruments made of aluminum (or other light metal) are not suitable for cleaning using a disinfectant RDG</li> </ul> </li> <li>› Please follow the manufacturer's instructions and the maintenance specifications for using of the autoclave and the disinfectant.</li> <li>› The method of use, reaction time and suitability of disinfectants and cleaning agents for certain types of instrument are covered by the manufacturers' instructions.</li> </ul>	
<p>These is a comprehensive discription of our detailed reprocessing advices. This can be found on our website, in our catalogue on page 218 and on request we'd be pleased to send you the printed information.</p>	

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


*Dica Diamonds basically consist of the shank plus a hardened, stainless profiled body which is coated with only selected, diamond grains using state-of-the-art electroplating techniques. This guarantees an extremely homogeneous, secure and long-lasting diamond coating. A choice of up to seven different grain sizes is also available for an even better result.*

*The use of coarse-grained diamonds (ISO 534, 544 and 554) can lead to increased heat generation. When using these products, special care should therefore be taken to ensure sufficient cooling and minimal application force. After using these diamond instruments, finishing should be performed in order to obtain an optimal surface roughness.*

*Instruments from ISO size 031 should be used with additional cooling (water spray).*

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
*Depth cuts may be used for crown, veneer or onlay preparation. Depth guides insure adequate but not too much tooth reduction.*

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

Flat end tapered or parallel diamonds with a rounded corner. Creates a modified shoulder with a preparation with a rounded internal angle and a 90° margin. Preferred margin for all-ceramic crowns.

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
This preparation is used for metal ceramic crowns or cast gold crowns. It creates a 45°-60° angle.

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Round end tapered or cylindrical diamond instruments create a chamfered preparation. Preferred preparation design for metal-ceramic crowns.

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*Provides ideal contour and shape for optimal interpretation by the scanning devices used in CAD/CAM fabrication or to take a precise impression.*

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
*Mainly used for design of stainless steel or cast gold crowns.*

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*An anatomically prepared occlusal surface provides adequate clearance and uniform thickness of the restoration. Football and wheel shape diamonds reproduce lingual concavity in crown preparation.*

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


*Refine shoulder preparation for improved crown adaption. Flatten the floor in inlay/onlay preparation.*

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

*Removing residual plaque from subgingival areas, Polishing the rootsurface, Removing subgingival calculus*

*Enlever la plaque subgingival, Polir la surface des racines, Enlever le calcul subgingival*

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
*Frequently used for crown preparation, to cut depth guides, occlusal rest preparation, endodontic access and cavity preparation.*

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
*Used for recreation of occlusal anatomy on restorations and rapid removal of obsolete restorations.*

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
*Special to improve conservative preparation design. Precise material removal due to small working part and slender neck. Minimally invasive shaping of cavities and maximum conservation of health tooth structure.*

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*Safe tip – Instruments with a non-cutting end*

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
*Instruments for embrasure contouring*

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
*Flat end diamonds with a rounded corner create rounded internal angles for ceramic or composite inlay and onlay preparation.*

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### 888L

Nadelform, kurz,  
langer Hals  
*needle-shaped, short,  
long neck*

FIG	SHANK	ISO	Ø
888LG	FG	806 314 539 534 ...	010
L mm			3,0
📏			5
APPLICATION			

### 889L

Nadelform,  
langer Hals  
*needle-shaped, long  
neck*

FIG	SHANK	ISO	Ø	
889L	FG	806 314 540 524 ...	009	010
889LG	FG	806 314 540 534 ...	009	010
889LF	FG	806 314 540 514 ...	009	010
L mm			3,5	4,0
📏			5	5
APPLICATION				

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
### MICRO DIAMONDS - MINIMALLY INVASIVE TREATMENT

MICRO DIAMANTEN - MINIMAL INVASIVE BEHANDLUNG | MICRO DIAMANTS - TRAITEMENT PEU INVASIF |  
MICRO DIAMANTES - TRATAMIENTO MÍNIMAMENTE INVASIVO

## 137

zylindrisch,  
Stirn rund  
*cylindrical, end  
hemispherical*

	FIG	SHANK	ISO	Ø
	137	FG	806 314 137 524 ...	007
L mm				2,1
				5
APPLICATION				

## 138

zylindrisch,  
Stirn rund  
*cylindrical, end  
hemispherical*

	FIG	SHANK	ISO	Ø
	138	FG	806 314 138 524 ...	007
L mm				3,6
				5
APPLICATION				

## 194

konisch, Stirn  
rund  
*conical, end  
hemispherical*

	FIG	SHANK	ISO	Ø
	194	FG	806 314 194 524 ...	007
L mm				2,1
				5
APPLICATION				

**271**

zylindrisch,  
Stirn rund  
cylindrical, end  
hemispherical

FIG	SHANK	ISO	Ø
271	FG	806 314 271 524 ...	007
L mm			2,1
			5
APPLICATION			

**277**

Ei  
egg

FIG	SHANK	ISO	Ø
277	FG	806 314 277 524 ...	009
L mm			1,1
			5
APPLICATION			

**295**

Torpedo, konisch  
torpedo, conical

FIG	SHANK	ISO	Ø
295	FG	806 314 295 524 ...	007
L mm			3,6
			5
APPLICATION			

**540**

Nadelform,  
langer Hals  
needle-shaped, long  
neck

FIG	SHANK	ISO	Ø
540	FG	806 314 540 524 ...	008
L mm			3,6
			5
APPLICATION			

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**697**

Rund, extra  
langer Hals  
spherical (round), extra  
long neck

FIG	SHANK	ISO	Ø
697	FG	806 314 697 524 ...	007
L mm			0,7
			5
APPLICATION			

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Corwn cutter for the separation of zirconia crowns and bridges

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
**CROWN CUTTERS - FOR HARD PRECIOUS METAL AND NON-PRECIOUS METAL ALLOYS**

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

**DIAMOND POLISHER FOR ZIRCONIUM OXIDE**

DIAMANTPOLIERER FÜR ZIRKONOXID | POLISSOIRS DIAMANTÉS POUR OXIDE DE ZIRCONIUM | PULIDOR DE DIAMANTE PARA ZIRCONIA

DIAMANTPOLIERER <i>Diamond polishers / Polissoirs diamantés / Pulidores diamantados</i>	DUNKELGRÜN <i>dark green / vert foncé / verde oscuro</i>	HELLGRÜN <i>light green / vert clair / verde claro</i>	ANWENDUNG & HYGIENE <i>Application &amp; Hygiene / Application &amp; Hygiène / Aplicación e higiene</i>	DREHZAHL <i>Rotation speed / Vitesse de rotation / Número de revoluciones</i>
Zweistufiges Diamantpoliersystem zum Abtragen, Glätten und Polieren von Zirkondioxid <i>two-step diamond polishing system for material removal, smoothing and gloss polishing of zirconium oxide / Système de polissoirs diamantés à deux étapes, pour l'enlèvement, le lissage et pour le lustrage de zirconium dioxyde / sistema de pulidores diamantados de dos pasos, para remover, alisar y pulir de circonio dióxido</i>	Glätten und Vorpulieren <i>smoothing and pre-polishing / lissage et pré-polissage / pulido inicial</i>	Hochglanzpolieren <i>gloss polishing / polissage brillant / pulido inal de alto brillo</i>		7'000 - 12'000 rpm max. 20'000 rpm

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
> *Polisher for removing, contouring, smoothing and high-lustre polishing of amalgam and gold.*  
> *Black (coarse): removal / contouring*  
> *Brown (medium): smoothing | Green (fine): high-lustre polishing*

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**POLISHERS FOR POLISHING MACROCOMPOSITES**

POLIERER ZUR POLITUR VON MAKROKOMPOSITEN | POLISSOIRS POUR LE POLISSAGE DE COMPOSITES MACRO | PULIDOR PARA EL PULIDO DE MACROCOMPOSITES

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


> Ceramic substitute for a scalpel used for exposing healed implants of adhesive residue

> Gingiva trimmer for gentle dilatation of the gingiva

> **Always used without any water-/spray-cooling!**

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

**1RF**

Fräser –  
Stahl rostfrei  
cutter - stainless steel

FIG	SHANK	ISO						
1RF	HP	330 104 001 001 ...	009	010	012	014	018	023
			5	5	5	5	5	5
APPLICATION								

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**Recommendation**

> Rotate the trephine a few times anti-clockwise to scribe a line. Then cut to the full depth intermittently in a clockwise direction.

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2550 - CROWN & BRIDGE SET

FOR CROWN & BRIDGE CHAMFER PREPARATION

		FIG	SHANK	ISO	Ø
		833	FG	806 314 277 524 ...	023
		833F	FG	806 314 277 514 ...	023
		850	FG	806 314 198 524 ...	018
		850F	FG	806 314 198 514 ...	018
		850	FG	806 314 198 524 ...	021
		850F	FG	806 314 198 514 ...	021
		852	FG	806 314 199 524 ...	014
		852F	FG	806 314 199 514 ...	014
		859L	FG	806 314 167 524 ...	010

2555 - PROFESSIONAL VENEER PREP


		FIG	SHANK	ISO	Ø
		833	FGL	806 315 277 524 ...	023
		833F	FGL	806 315 277 514 ...	023
		815	FG	806 314 042 524 ...	022
		815	FG	806 314 042 524 ...	026
		852	FG	806 314 199 524 ...	012
		852F	FG	806 314 199 514 ...	012
		852	FG	806 314 199 524 ...	016
		852F	FG	806 314 199 514 ...	016
		890EF	FG	806 314 160 504 ...	010

**ZIRCONIUM OXIDE PREPARATION**

ZIRKONDIOXID BEARBEITUNG | PRÉPARATION DU OXYDE DE ZIRCONIUM | PREPARACIÓN DEL ÓXIDO DE CIRCONIO

**2512 - ZIRCONIUM OXIDE  
HIGH GLOSS POLISHING**

		FIG	SHANK	ISO	∅
		ZIR9861M	RA	803 204 243 521 ...	040
		ZIR9861F	RA	803 204 243 511 ...	040
		ZIR9862M	RA	803 204 030 521 ...	060
		ZIR9862F	RA	803 204 030 511 ...	060
		ZIR9863M	RA	803 204 543 521...	140
		ZIR9863F	RA	803 204 543 511 ...	140

**2510 - ZIRCONIUM OXIDE  
PREPARATION & FINISHING**

		FIG	SHANK	ISO	∅
		Z833	FG	806 314 277 314 ...	023
		Z850	FG	806 314 198 314 ...	018
		GP662	RA	655 204 288 526 ...	035
		ZIR9861M	RA	803 204 243 521 ...	040
		ZIR9861F	RA	803 204 243 511 ...	040

ALL CERAMIC PREPARATION

VOLLKERAMIK BEARBEITUNG | TOUT CÉRAMIQUE PRÉPARATION | PREPARACIÓN DEL TODO CERÁMICA

2522 - ALL CERAMIC  
HIGH GLOSS POLISHING

		FIG	SHANK	ISO	Ø
		LS9871M	RA	803 204 243 522 ...	040
		LS9871F	RA	803 204 243 512 ...	040
		LS9872M	RA	803 204 030 522 ...	060
		LS9872F	RA	803 204 030 512 ...	060
		LS9873M	RA	803 204 543 522 ...	140
		LS9873F	RA	803 204 543 512 ...	140


2520 - ALL CERAMIC  
PREPARATION & FINISHING

		FIG	SHANK	ISO	Ø
		833F	FG	806 314 277 514 ...	023
		850	FG	806 314 198 514 ...	018
		SD662F	RA	805 204 297 514 ...	030
		LS9871M	RA	803 204 243 522 ...	040
		LS9871F	RA	803 204 243 512 ...	040

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*Specially developed for preparing zirconia using a laboratory turbine and water cooling*

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- > Before using for the first time the instruments must be cleaned on a cleaning stone, such as the 529D, to remove any production residue. The cleaning stone should be moistened beforehand.
- > Sinter diamonds should be cleaned periodically using the cleaning stone.
- > Avoid applying excessive pressure to sinter diamonds, as they break much more easily than electroplated diamonds..

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

**DIAMOND DISCS "FLEXIBLE" - LARGE DIAMETER**

DIAMANTSCHIBEN "FLEXIBEL" GROSSER DURCHMESSER | DISQUES DIAMANTÉS „FLEXIBLES“ GRAND DIAMÈTRE | DISCOS DE DIAMANTE "FLEXIBLES" GRAN DIÁMETRO

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**DIAMOND DISCS "SUPERFLEXIBLE" ONE-SIDED**

DIAMANTSCHIBEN "SUPERFLEXIBEL" EINSEITIG | DISQUES DIAMANTÉS „SUPER FLEXIBLES“ D'UN CÔTÉ | DISCOS DE DIAMANTE "SUPERFLEXIBLES" UN LADO

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




> This cylindrical cutter is particularly thin and has a special cutting geometry. It is ideal for trimming restorations from CAD/CAM blanks. It glides through the material without tearing or build-up of pressure. The extra long shank provides optimal access, even with high blank thicknesses.

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> for preparing acrylic and vacuum-formed templates

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> for efficient preparation of all hard materials  
> for producing very smooth surfaces

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> for rapid and effective trimming of acrylics and also plaster

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> for the finest preparation of metal alloys and composites

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> mainly for finer preparation of precious metals, non-precious metals, acrylics and porcelain

> moderate, highly selective material removal, smooth workpiece surface

> instruments with an M in the Match Code have a shorter working section to allow more detailed preparation

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> for preparing larger surfaces of precious metals, non-precious metals and acrylics  
 > efficient material removal, smooth workpiece surfaces, smooth running of the tool  
 > very universal application options, consequently instruments are changed less frequently  
 > instruments with an M in the Match Code have a shorter working section to allow more detailed preparation

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> for rough, efficient pre-preparation of large surfaces of precious metals, non-precious metals and acrylics  
> for finer trimming of plaster

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> ideal for trimming plaster, also suitable for very rough preparation of acrylics  
> very high material removal  
> large chip escape channels prevent clogging with grinding chips

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> better protection of the blade against the edges breaking off, which considerably increases the service life of the cutters  
> the cutters also have perceptibly better concentricity without chatter  
> for preparing all hard metal alloys such as chrome-cobalt, non-precious metals and semi-precious alloys

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> for effective trimming of permanently soft relining material

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> special cut for quick, effective preparation of titanium and titanium alloys

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
>with this specially designed fissure instrument, the fissures in the occlusion of pre-sintered zirconia restorations can be finished to a high esthetic level.

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**For the highly esthetic material zircon oxide simply, reliably and precisely**

With Dendia Zirojet, reconstruction using zircon oxide, this increasingly valued high-tech material, is unbelievably easy and precise. The outstanding properties of zircon oxide – extreme wear-resistance, smooth surface and minimal heat conductivity – fulfil at the same time the toughest and most sensitive demands for high end finishing.

**Fast and precise**

- Sensitive, reliable and extremely precise grinding
- Simpler handling, better levelling, less contact pressure
- No formation of microcracks during dry processing
- Ideal for dealing with edges-
- Little heat development, no sparking
- Low revolutions per minute – high wear- resistance
- Central, vibration-free movement

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*The 100% diamond abrasive opens up completely new perspectives under well-established operating conditions: very fast, precise handling and thus high performance and rapid workflow. And all this with a previously unknown level of wear-resistance and thus higher profitability.*

**High service life, universal application**

- Optimised solution for the use of ceramics, pressed ceramics and synthetic materials
- Very good cutting results, particularly for glass ceramics and synthetic materials
- Low heat development
- Unbeatable service life
- Universal application
- Highest wear-resistance
- Reliable workmanship

FARBE <i>Colour   couleur   colores</i>		KÖRNUNG <i>Grain   Grain   Grano</i>	
	roter Ring <i>red ring   bague rouge   anillo rojo</i>	fein <i>fine   fin   fino</i>	= F

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
*Silicon carbide for preparing porcelain and acrylics*

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




*High-Grade corundum for preparing medium-hard metal alloys and chrome-cobalt*

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


*High-Grade corundum for preparing semi-precious and non-precious metal alloys*

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

*High-Grade corundum for  
preparing precious metal alloys*

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
*High-Grade corundum for preparing hard metal alloys and non-precious metal alloys*

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*The bend-resistant nylon cut-off wheels with low abrasion, consisting of a phenolic resin base with standard corundum abrasive, are ideal for cutting off and preparing chrome-cobalt alloys (CrCo) as well as precious and non-precious metal alloys.*

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*The nylon separating discs, consisting of a phenolic resin and silicon carbide, are ideal for separating and preparing porcelain with minimal heat build up and without contaminating grinding marks.*

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
SILICONE POLISHERS FOR PREPARING DENTURE ACRYLIC

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
Flexible polisher for smoothing precious and non-precious metal alloys

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**DIAMOND POLISHER FOR COMPOSITES & NANO- / HYBRID CERAMICS**

DIAMANTPOLIERER ZUR KOMPOSIT & NANO- / HYBRIDKERAMIK BEARBEITUNG | POLISSOIRS DIAMANTÉS POUR L'USINAGE DE COMPOSITES & NANO CÉRAMIQUE | PULIDORES DIAMANTADOS PARA TRABAJOS DE COMPOSITES & NANO CERÁMICA

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*One-step diamond polisher for smoothing and polishing of composite and nano- / hybrid-ceramic*

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DIAMOND POLISHER FOR ZIRCONIUM OXIDE

DIAMANTPOLIERER FÜR ZIRKONOXID | POLISSOIRS DIAMANTÉS POUR OXIDE DE ZIRCONIUM | PULIDOR DE DIAMANTE PARA ZIRCONIA

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**Diamond polishers - ZIR Gloss**

> two-step diamond polishing system for material removal, smoothing and polishing of zirconium oxide  
 > dark green: smoothing and pre-polishing / light green: gloss polishing

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

**Diamond polishers - LS Gloss**

> two-step diamond polishing system for material removal, smoothing and polishing of lithium disilicate and feldspar ceramic  
 > red-brown: smoothing and pre-polishing / beige: gloss polishing

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
*after use on metals, clean the polishers to prevent striation marks forming on the porcelain*

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

*Soft goat hair brush - for polishing denture acrylics, composites & precious metal alloys with polishing paste*

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*Horsehair medium - for polishing non precious metal crown and bridge restorations*

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*White Chungking bristles hard - for polishing of composite materials white polishing paste*

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*Non-woven fabric with medium grit/ double layered - for preparing the surfaces of soft relines, vacuum-formed templates, acrylic dentures and custom impression trays*

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*Non-woven fabric with fine grit/ double layered - for inishing the surfaces of soft relines, vacuum-formed templates, acrylic dentures and custom impression trays*

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*Miniature brush made from fine cotton yarn - shine polishing of plastic prosthetics, composites, crown and bridge work made of precious metal, NEM and CoCr with emulsion*

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*Miniature buffing wheel made of fine cloth - high-luster polishing of acrylic dentures, composites, crown and bridge work made of precious metal, non-precious metal and CoCr with polishing paste*

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*Miniature buffing wheel made of microfibre leather - shine polishing of plastic prosthetics, composites, crown and bridge work made of precious metal, NEM and CoCr with polishing paste*

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*Abrasive-impregnated brushes for polishing denture acrylics, composites and precious metal alloys*

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
2410 - ZIRCONIUM OXIDE  
PREPARATION PROFESSIONAL


		FIG	SHANK	ISO	∅
		C31RL	HPL	500 105 143 007...	010
		9310M	HP	653 104 292 524...	050
		863	HP	806 104 250 524...	012
		515L	HP	310 104 467 211...	012
		SZ652R	HP	655 104 199 526...	035
		ZIR9868M	HP	803 104 543 521...	140
		ZIR9868F	HP	803 104 543 511 ...	140
		ZIR9865M	HP	803 104 303 521...	260
		ZIR9865F	HP	803 104 303 511 ...	260

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## INSTRUCTIONS FOR THE PROCESSING (CLEANING, DISINFECTION, AND STERILIZATION) OF INSTRUMENTS from DENDIA GmbH Issued: August 2017

The medical devices produced and sold by DENDIA GmbH are re-usable unless their label contains explicit information to the contrary. However, as a rule, it is the sole responsibility of the doctor/expert using the devices to decide whether, depending on the respective case and the potential wear and tear of the products, he can re-use the products and how frequently he uses them. In case of doubt, it is always advisable to discard the products early and to replace them. The manufacturer DENDIA GmbH cannot guarantee the faultless function and performance of the products combined with a maximum degree of safety if the products are overused. These reprocessing instructions apply in principle to all medical devices making up the product range supplied by DENDIA GmbH. Any particular features and/or exclusions that only concern individual items or groups of items are referred to separately.

### FUNDAMENTAL POINTS

All instruments are to be cleaned, disinfected, and sterilized prior to each application; this is required as well for the first use after delivery of the unsterile instruments (cleaning and disinfection after removal of the protective packaging, sterilization after packaging). An effective cleaning and disinfection is an indispensable requirement for an effective sterilization of the instruments.

You are responsible for the sterility of the instruments. Therefore, please ensure that only sufficiently device and product specifically validated procedures will be used for cleaning, disinfection, and sterilization, that the used devices (WD, sterilizer) will be maintained and checked regularly, as well as that the validated parameters will be applied for each cycle.

Please pay attention to avoid a higher contamination of the complete bur block during application; otherwise it is necessary to clean and disinfect the bur block as well as all instruments inside (after removal).

Additionally, please pay attention to the legal provisions valid for your country as well as to the hygienic instructions of the doctor's practice or of the hospital. This applies particularly to the different guidelines regarding the inactivation of prions (not relevant for USA).

**SOME INSTRUMENTS REQUIRE ADDITIONAL ASPECTS. FOR THIS, PAY ATTENTION TO CHAPTER "SPECIFIC ASPECTS"**

### CLEANING AND DISINFECTING

#### BASIC:

If possible, an automated procedure (WD (Washer-Disinfector)) should be used for cleaning and disinfection of the instruments. A manual procedure – even in case of application of an ultrasonic bath – should only be used if an automated procedure is not available; in this case, the significantly lower efficiency and reproducibility of a manual procedure has to be considered.

The pre-treatment step is to be performed in both cases.

#### PRE-TREATMENT:

Please remove coarse impurities of the instruments directly after application (within a maximum of 2 h).

Procedure:

1. Rinse the instruments at least 1 min under running water (temperature < 35 °C/95 °F).
2. Soak the instruments at least for the given soaking time in the pre-cleaning solution<sup>1</sup> (by the use of an ultrasonic bath) so that the instruments are sufficiently covered. Pay attention that there is no contact between the instruments. Assist cleaning by careful brushing with a soft brush (at least three times after at beginning of soaking, aids see chapter "Specific aspects").
3. Activate ultrasonic treatment for an additional soaking time (but not less than 5 min).
4. Then, remove the instruments of the pre-cleaning solution and post-rinse them at least three times intensively (at least 1 min) with water.
5. In case of still visible contamination repeat steps 2, 3, and 4, otherwise discard the instrument. This is especially relevant for diamond instruments.

Pay attention to following points during selection of the cleaning detergent:

- > fundamental suitability for the cleaning of instruments made of metallic or plastic material
- > suitability of the cleaning detergent for ultrasonic cleaning (no foam development)
- > compatibility of the cleaning detergent with the instruments (see chapter „material resistance.“)

Pay attention to the instructions of the detergent manufacturer regarding concentration, temperature and soaking time as well as post-rinsing. Please use only freshly prepared solutions as well as only demineralized sterile or low contaminated water (max. 10 germs/ml) as well as low endotoxin contaminated water (max. 0.25 endotoxin units/ml), for example purified/highly purified water, and a soft, clean, and lint-free cloth and/or filtered air for drying, respectively.

<sup>1</sup> In case of application of a cleaning and disinfection detergent for this (e.g. in consequence of personnel's safety) please consider, that this should be aldehyde-free (otherwise fixation of blood impurities), possess a fundamentally approved efficiency (for example VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking), be suitable for the disinfection of instruments made of metallic or plastic material, and be compatible with the instruments (see chapter „material resistance.“). Please consider, that a disinfectant used in the pre-treatment step serves only the personnel's safety, but cannot replace the disinfection step later to be performed after cleaning.

### AUTOMATED CLEANING/DISINFECTION (WD (WASHER-DISINFECTOR)):

Pay attention to following points during selection of the WD:

- > fundamentally approved efficiency of the WD (for example CE marking according to EN ISO 15883 or DGHM or FDA approval/clearance/registration)
- > possibility for an approved program for thermal disinfection (A0 value > 3000 or – in case of older devices – at least 5 min at 90°C/194 °F; in case of chemical disinfection danger of remnants of the disinfectant on the instruments)
- > fundamental suitability of the program for instruments as well as sufficient rinsing steps in the program
- > post-rinsing only with demineralized sterile or low contaminated water (max. 10 germs/ml, max. 0.25 endotoxin units/ml), for example purified/highly purified water
- > only use of filtered air (oil-free, low contamination with microorganisms and particles) for drying
- > regularly maintenance and check/calibration of the WD

Pay attention to following points during selection of the cleaning detergent:

- > fundamental suitability for the cleaning of instruments made of metallic or plastic material
- > additional application – in case of non-application of a thermal disinfection – of a suitable disinfectant with approved efficiency (for example VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking) compatible to the used cleaning detergent
- > compatibility of the used detergents with the instruments (see chapter „material resistance.“)

Pay attention to the instructions of the detergent manufacturers regarding concentration, temperature and soaking time as well as post-rinsing.

#### PROCEDURE:

1. Transfer the instruments in the WD by the use of a small pieces basket.
2. Start the program.
3. Remove the instruments of the WD after end of the program.
4. Check and pack the instruments immediately after the removal (see chapters „check.“, „maintenance.“, and „packaging.“, if necessary after additional post-drying at a clean place).

The fundamental suitability of the instruments for an effective automated cleaning and disinfection was demonstrated by an independent, governmentally accredited and recognized [§ 15 (5) MPG] test laboratory by application of the WD G 7836 CD, Miele & Cie. GmbH & Co., Gütersloh, (thermal disinfection) and the pre-cleaning and cleaning detergent Neodisher mediclean forte (Dr. Weigert GmbH & Co. KG, Hamburg) considering to the specified procedure.

#### MANUAL CLEANING AND DISINFECTION:

Pay attention to following points during selection of the cleaning and disinfection detergents:

- > fundamental suitability for the cleaning and disinfection of instruments made of metallic or plastic material
- > in case of application of an ultrasonic bath: suitability of the cleaning detergent for ultrasonic cleaning (no foam development)
- > application of a disinfectant with approved efficiency (for example VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking) compatible with the used cleaning detergent
- > compatibility of the used detergents with the instruments (see chapter „material resistance.“)

Combined cleaning/disinfection detergents should not be used.

Only in case of extremely low contamination (no visible impurities) combined cleaning/disinfection could be used.

Pay attention to the instructions of the detergent manufacturers regarding concentration, temperature and soaking time as well as post-rinsing. Please use only freshly prepared solutions as well as only demineralized sterile or low contaminated water (max. 10 germs/ml) as well as low endotoxin contaminated water (max. 0.25 endotoxin units/ml), for example purified/highly purified water, and a soft, clean, and lint-free cloth and/or filtered air for drying, respectively.

#### PROCEDURE:

##### CLEANING

1. Soak the instruments for the given soaking time in the cleaning solution (by the use of an ultrasonic bath) so that the instruments are sufficiently covered. Pay attention that there is no contact between the instruments. Assist cleaning by careful brushing with a soft brush (at least three times after at beginning of soaking, aids see chapter "Specific aspects").
2. Activate ultrasonic treatment for an additional soaking time (but not less than 15 min).
3. Then, remove the instruments of the cleaning solution and post-rinse them at least three times intensively (at least 1 min) with water.
4. Check the instruments (see chapters „check.“ and „maintenance.“).

##### DISINFECTION

5. Soak the instruments for the given soaking time in the disinfectant solution so that the instruments are sufficiently covered. Pay attention that there is no contact between the instruments.
6. Then, remove the instruments of the disinfectant solution and post-rinse them at least five times intensively (at least 1 min) with water.
7. Dry and pack the instruments immediately after the removal (see chapter „packaging.“, if necessary after additional post-drying at a clean place).

The fundamental suitability of the instruments for an effective cleaning and disinfection was demonstrated by an independent, governmentally accredited and recognized [§ 15 (5) MPG] test laboratory by application of the pre-cleaning and cleaning detergent Cidezme/Enzol and the disinfectant Cidex OPA (Johnson & Johnson GmbH, Norderstedt) considering the specified procedure.

#### CHECK

Check all instruments after cleaning or cleaning/disinfection, respectively, on corrosion, damaged surfaces, and impurities. Do not further use damaged instruments (for limitation of the numbers of re-use cycles see chapter „reusability.“). Still dirty instruments are to be cleaned and disinfected again.

#### MAINTENANCE

Instrument oils or grease must not be used with the exception of steel instruments. In that case use only instrument oils (white oil) admitted to steam sterilization considering the maximum possible sterilization temperature, with approved biocompatibility and without mono-, di, or triethanolamine as corrosion inhibitor.

#### PACKAGING

Please insert the cleaned and disinfected instruments in the corresponding bur blocks (if required) and pack them in single-use sterilization packagings (single or double packaging), which fulfill the following requirements (material/process):

- > EN ISO/ANSI AAMI ISO 11607 (for USA: FDA clearance)
- > suitable for steam sterilization (temperature resistance up to at least 142 °C [288 °F], sufficient steam permeability)
- > sufficient protection of the instruments as well as of the sterilization packagings to mechanical damage

#### STERILIZATION

Please use for sterilization only the listed sterilization procedures; other sterilization procedures must not be applied.

##### Steam sterilization

- > fractionated vacuum/dynamic air removal procedure<sup>2,3</sup> (with sufficient product drying<sup>5</sup>)
- > steam sterilizer according to EN 13060/EN 285 or ANSI AAMI ST79 (for USA: FDA clearance)
- > validated according to EN ISO 17665 (valid IQ/OQ (commissioning) and product specific performance qualification (PQ))
- > maximum sterilization temperature 138 °C (280 °F; plus tolerance according to EN ISO 17665)
- > sterilization time (exposure time at the sterilization temperature):

AREA	FRACTIONATED VACUUM/DYNAMIC AIR REMOVAL	GRAVITY DISPLACEMENT
USA	at least 4 min at 132 °C (270 °F), drying time at least 20 min <sup>4</sup>	not recommended
other countries	at least 3 min <sup>4</sup> at 132 °C (270 °F) / 134 °C (273 °F), drying time at least 20 min <sup>4</sup>	not recommended

<sup>2</sup> at least three vacuum steps

<sup>3</sup> The less effective gravity displacement procedure must not be used in case of availability of the fractionated vacuum procedure, will require significantly longer sterilization times and is to be validated dependent on product, packaging, sterilizer, program, and parameters under sole responsibility of the user.

<sup>4</sup> The effectively required drying time depends directly on parameters in sole responsibility of the user (load configuration and density, sterilizer conditions, ...) and by this is to be determined by the user. Nevertheless, drying times less than 20 min must not be applied.

<sup>5</sup> respectively 18 min (inactivation of prions, not relevant for USA)

The fundamental suitability of the instruments for an effective steam sterilization was demonstrated by an independent, governmentally accredited and recognized [§ 15 (5) MPG] test laboratory by application of the steam sterilizer HST 6x6x6 (Zirbus technology GmbH, Bad Grund) and the fractionated vacuum/dynamic air removal procedure. For this, typical conditions in clinic and doctor's practice as well as the specified procedure were considered.

The flash/immediate use sterilization procedure must not be used. Do not use dry heat sterilization, radiation sterilization, formaldehyde and ethylene oxide sterilization, as well as plasma sterilization.

#### STORAGE

Please store the instruments after sterilization in the sterilization packagings at a dry and dust-free place.

#### MATERIAL RESISTANCE

Please take care that the listed substances are not ingredients of the cleaning or disinfection detergent:

- > organic, mineral, and oxidizing acids (minimum admitted pH-value 5.5)
- > strong lyes (maximum admitted pH-value 11, neutral/enzymatic or alkaline cleaner recommended<sup>6</sup>)
- > organic solvents (for example: acetone, ether, alcohol, benzene)
- > oxidizing agents (for example: hydrogen peroxide)
- > halogens (chlorine, iodine, bromine)
- > aromatic, halogenated hydrocarbons

<sup>6</sup> For the bur blocks alkaline cleaners must not be applied (maximum admitted pH-value 9).

Please do not clean any instruments and bur blocks by use of metal brushes or steel wool.

Please do not expose any instruments and bur blocks to temperatures higher than 142 °C (288 °F)!

Please do not apply acidic neutralizing agents or cleaning acids.

#### REUSABILITY

The instruments can be reused – in case of adequate care and if they are undamaged and clean as indicated in chapter "Specific aspects". The user is responsible for each further use as well as for the use of damaged and dirty instruments (no liability in case of disregard).

# APPLICATION AND SAFETY INSTRUCTIONS

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## ATTACHMENT A: SPECIFIC ASPECTS

Diamond products and ceramic grinding tools:

- > Use particular care during the cleaning of the grinding surfaces and ensure that all residues are removed

Bur blocks/instrument trays:

- > Cleaning and disinfecting only without products being loaded (products must be cleaned and disinfected whilst they are in the bur block/ instrument tray)

i
Das Dokument in der jeweils aktuellsten und gültigen Version finden sie auf unserer Webseite <a href="http://www.dendidental.com">www.dendidental.com</a>
The document in the latest and applicable version you find on the website <a href="http://www.dendidental.com">www.dendidental.com</a>
Vous trouverez les documents dans la dernière version sur notre site <a href="http://www.dendidental.com">www.dendidental.com</a>
Documentos en la versión más reciente encontrase en nuestra página web <a href="http://www.dendidental.com">www.dendidental.com</a>

Instrument group	brush	specific/additional procedure in case of						maximum admitted cycle number (confirmed by validation, but dependent on specific application)	recommended classification according to KRINKO/RKI/BfArM guidance (only German, with respect to intended use)
		pretreatment	manual cleaning/ disinfection	automated cleaning/ disinfection	maintenance	packing	sterilization		
stainless steel instruments	standard	standard	standard	standard	lubrication <u>not</u> admitted	standard	standard	10	critical B
regular steel instruments	standard	standard	standard	standard	lubrication <u>recommended</u>	standard	standard	10	critical B
silicone polisher	standard	standard	standard	standard	lubrication <u>not</u> admitted	standard	standard	5	critical B
endodontic instruments without stopper	endodontic brush	standard	standard	standard	lubrication <u>not</u> admitted	use of bur blocks not admitted	use of bur blocks not admitted	10	critical B
endodontic instruments with stopper	endodontic brush	mounted	mounted move the stopper at least three times during disinfection	mounted	lubrication <u>not</u> admitted	use of bur blocks not admitted	use of bur blocks not admitted	10	critical B
all other instruments	standard	standard	standard	standard	lubrication <u>not</u> admitted	standard	standard	10	critical B

## GENERAL APPLICATION AND SAFETY INSTRUCTIONS

for the medical device from DENDIA GmbH  
Issued: August 2017

- > DENDIA GmbH products (dental, maxillary surgery, general surgery,) must only be used by dentists, doctors and/or the respective experts who, due to their training and experience, are intensely familiar with the use of these products and who have the corresponding expertise in the respective specialist fields. The use of surgical products requires relevant expertise and experience in dental implantology, maxillary surgery and/or other surgical fields including diagnosis, preoperative planning and surgical techniques.
- > It is the sole responsibility of the doctor in charge who, depending on the respective situation (indication), decides on the actual use of the DENDIA GmbH products for each individual case
- > All DENDIA GMBH products have been developed for specific applications. Therefore, inappropriate use can result in the premature wear and tear of the products and put patients and users at risk.

## APPLICATION

- > In order to avoid damaging the instruments, they must be removed from the blister pack by pulling off the back-sheet.
- > It is essential to only use turbines as well as hand and angle pieces that are technically and hygienically faultless, maintained and cleaned.
- > The instruments must be rotating when applied on material. They should not be placed on material and then brought to rotation.
- > Rotating instruments need to be clamped as far down as possible with their speed set before applying them on the object. are used with the rotary instruments.
- > Using the instruments for canting or leveraging should be avoided as it increases the risk of breakage.
- > Depending on the application, it is recommended to use protective goggles while using the instruments. Users of diamond disks should use a disk protection device.
- > Inappropriate use of the products leads to badly executed work and increased risk.
- > When working with dry materials, it is recommended to use a suction cleaning device.
- > In particular, users of hand tools should take care to use them gently and with consideration.
- > The user must at all times avoid touching the instruments and parts without protection (protective gloves should be worn).
- > Thermal bone damage caused by rotating and oscillating tools (e.g. pilot burr, conical burr, expansion burr) must at all times be avoided (user training, working at low speed and with sufficient cooling).
- > During intraoral application attention has to be made to the fact that the products are protected against aspiration or falling on the floor.

## USE OF PRESSURE

- > Users of the instruments should at all times avoid applying excessive pressure. This can damage the working part of the instruments and cause the cutting edges to break off. At the same time, it generates excessive heat.
- > The use of excessive pressure when using grinding tools can cause the abrasive particles to break off or the instrument to become clogged and lead to heat generation.
- > During polishing, excess pressure can lead to heat generation.
- > Due to overheating, excess pressure can damage the dental pulp or, due to broken off cutting edges, it can result in undesired rough surfaces. In such cases, even instrument breakage cannot be excluded.

## COOLING

- > In order to avoid excessive heat generation during preparation, a sterile water/sodium chloride solution supplied via a permanent external feeding device should be used to ensure sufficient cooling during use of the instruments.
- > When using FG instruments that are more than 22 mm long or whose head diameter exceeds 2 mm, additional external cooling is required.
- > Insufficient cooling will lead to irreversible damage to the bone and/or the adjacent tissue.

## STORAGE, DISINFECTION, CLEANING AND STERILIZATION

- > Unless there is explicit information to the contrary, all DENDIA GmbH products are supplied in non-sterile packaging and, depending on the application, they need to be sterilised prior to use. Prior to their first use on the patient and immediately after each use, all products need to be disinfected and sterilised. Inappropriate cleaning and sterilising of the instruments can result in the patient being infected with harmful bacteria.
- > You will find detailed instructions for the disinfecting, cleaning and sterilising of products in the Instructions for the processing of instruments produced by DENDIA GmbH on the previous page. We would also be happy to provide you with these instructions at your request. They are also available on the internet at [www.dendidental.com](http://www.dendidental.com).
- > The products should be stored in appropriate, hygienically maintained containers. The same applies to sterilised instruments. The stored products must be protected from dust, humidity and recontamination. Instructions as to maximum storage duration must be adhered to.

## SPEED RECOMMENDATIONS FOR ROTARY INSTRUMENTS

- > Following the instrument-specific speed recommendations produces the best results.
- > Exceeding the maximum admissible speed (rpm) when using long and pointed instruments tends to produce vibrations that can lead to the destruction of the instrument.
- > When using working parts with diameters exceeding the thickness of the shaft, excessive speed can release great centrifugal forces that may cause the shaft to bend and/or the instrument to break. Therefore, the maximum admissible rpm must never be exceeded. Please consult the manufacturer's information (see catalogue or [www.dendidental.com](http://www.dendidental.com)) for the recommended and the maximum admissible speed ranges. Non-compliance with the maximum

- > Generally, the following rules apply:
  - The larger the working part of an instrument the lower the speed
  - Surgical instruments: suitable for geared down micro-motor hand and angle pieces 10:1 with stable ball bearings. Speed 600 to 800 rpm with physical and, possibly, sterile external cooling or internal cooling when using the respective hand piece.

## DISCARDING WORN INSTRUMENTS AND PARTS

- > DENDIA GmbH products can principally be reused several times – unless specifically indicated and labelled otherwise. Rotating instruments are subject to wear. The option of and accountability for multiple use of a product and the frequency of application is solely the decision and own responsibility of the treating clinician based on the application in each case and the possible wear of the products. If in doubt, the products should always be sorted out early and replaced.
- > Broken off cutting edges of instruments cause vibrations and great forces of pressure, which, in turn, leads to broken preparation corners and rough surfaces.
- > Bare patches on diamond instruments indicate a lack of abrasive particles and can be a sign of blunt instruments. This leads to excessive temperatures during instrument use.
- > Instruments that are bent and/or do not run true should be discarded forthwith.
- > With the reuse of disposable products the risk of infection cannot be excluded and a risk-free functional safety cannot be guaranteed.

## ADDITIONAL INSTRUCTIONS REGARDING THE USE OF TREPANS

- > When using trepans, you have to proceed with particular care. For example, it is advisable not to exceed the recommended rpm speed ranges.
- > In order to prepare for the actual use of a trepan, it should be set to produce counter-clockwise rotations creating a groove in the bone. Afterwards the trepan can be inserted into this groove and, using clockwise rotations, it can be moved further down.
- > Carrying out a prior X-Ray is essential to establish the maximum possible drilling depth and to maintain the necessary distance, for example, to the mandibular nerve. As an additional safety measure to spare the nerve, the axial direction of the trepan countersink attachment, based on the sagittal level of the ascending branch, must be milled laterally at an angle of approx. 15-20°.

$$\varnothing < 045 \quad \begin{matrix} \curvearrowright \\ < 4000 \\ \text{min}^{-1} \end{matrix} \quad \varnothing \geq 045 \quad \begin{matrix} \curvearrowright \\ < 2000 \\ \text{min}^{-1} \end{matrix}$$

## FURTHER COMMENTS:

- > Due to statutory regulations, returned goods can, on principle, only be accepted if the complete batch number is provided. This number can be found on the product packaging.

# GENERAL DELIVERY AND PAYMENT TERMS AND CONDITIONS

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## 1. General / scope

- 1.1 Our delivery and payment terms and conditions are exclusive; terms to the contrary or deviating from our delivery and payment terms and conditions and those of the purchaser's terms and conditions which are not dealt with in the present delivery and payment terms and conditions will not be accepted by us unless we have given our express consent to their validity in writing. Our delivery and payment terms and conditions shall be applicable also in those cases where we effect delivery to the purchaser without reservation in the knowledge of the purchaser's terms and conditions being contrary to or deviating from our delivery and payment terms and conditions or not dealt with in our delivery and payment terms and conditions.
- 1.2 All agreements made between us and the purchaser for purposes of executing a contract are laid down in writing in this contract and in these simultaneously effective delivery and payment terms and conditions.
- 1.3 Our delivery and payment terms and conditions are applicable also to repeat orders.

## 2. Quotation / quotation documents

- 2.1 Our quotation is non-binding unless an agreement to the contrary has been made.
- 2.2 Illustrations, drawings, calculations and other product, application and project related documents shall remain our property and are subject to our copyright even if we pass them on to the purchaser; they must not be copied nor made available to third parties without our express prior consent in writing.

## 3. Prices / payment terms / minimum quantities

- 3.1 Insofar as no agreements to the contrary have been made, our prices are "ex works" exclusive of packaging and transport costs; these will be charged for separately.
- 3.2 Our prices do not include statutory VAT; this will be shown separately on the invoice at the statutory rate applicable at the date of invoice.
- 3.3 Minimum order quantities for international customers are as follows: wholesaler and export customers EUR 250.- order volume. The minimum order per each order position is EUR 20.- for regular catalogue articles. Special articles not shown in catalogue will be offered in reasonable lot sizes.
- 3.4 The deduction of discount requires special agreement in writing.
- 3.5 Insofar as no agreement to the contrary is made, the purchase price is payable net (without deductions) within 30 days from the date of invoice. Should there be a payment delay on the part of the purchaser we shall be entitled to charge interest for delay of 8% above the base rate published at the time.
- 3.6 The purchaser shall have the right to off-set payments only if its counter-claims have been recognised by declaratory judgment, are undisputed or have been acknowledged by us. The purchaser is entitled to exercise a right of retention insofar as its counter-claim is based on the same contractual relationship.

## 4. Delivery / delivery period

- 4.1 Commencement of the delivery period indicated by us presupposes the clarification of all technical questions.
- 4.2 Compliance with our delivery obligation moreover presupposes the timely and proper fulfilment of the purchaser's obligations, in particular compliance with the agreed payment terms. The plea of non-fulfilment of the contract remains subject to reservation.
- 4.3 Should the purchaser be in default of acceptance or infringe other participation obligations, we shall be entitled to require compensation for the loss, including possible additional expenditure, incurred. We reserve the right to make further claims.
- 4.4 Insofar as the preconditions of section 4.3 have been fulfilled, the risk of accidental loss or accidental perishing of the purchased item passes to the purchaser at the moment when the purchaser is in default of acceptance or debtor's delay.
- 4.5 Should non-compliance with the delivery period be attributable to force majeure, e.g. mobilisation, war, terror, riot or similar events, e.g. strike, lock-out, the delivery period shall be extended by the duration of the events causing the delay, insofar as the obstacles are demonstrably of not inconsiderable influence on the completion or delivery. This applies also if these circumstances arise for sub-contractors.
- 4.6 We shall be freed from our delivery obligation if we fail to be supplied ourselves in due course with the right goods ordered for the fulfilment of the contract, this not being our fault.
- 4.7 Partial deliveries are permissible within reasonable limits.
- 4.8 We are liable in the case of a delay of delivery in accordance with the statutory provisions.
- 4.9 We are also liable in accordance with the statutory provisions, insofar as the delay of delivery is due to a wilful or grossly negligent infringement of the contract attributable to us; any fault committed by our representatives or person employed in the performance of our obligations is attributable to us.
- 4.10 Insofar as the delay of delivery is due merely to simple negligence and liability is not mandatory as a result of loss of life, bodily injury or injury to health, our liability for losses due to delay is limited to the purchaser being able to demand, for each complete week of the delay, 0.5 percent, though in total a maximum of 5 percent, of the price of that part of the delivery which cannot be appropriately taken into operation as a result of the delay. This does not entail any modification of the burden of proof to the detriment of the purchaser. The purchaser's statutory right of rescission remains unaffected.
- 4.11 Should the delivery be delayed for reasons for which the purchaser is responsible, storage costs amounting to 0.5 percent of the price of the items supplied, though not more than 5 percent in total, can be charged to the purchaser for every month that has commenced. On completion of an appropriate additional period, we shall be entitled to dispose of the item supplied elsewhere and to supply the purchaser after an extended period. The contractual parties are at liberty to provide proof that higher, lower or no storage costs have arisen. The statutory rights to withdraw from the contract and to demand compensation remain unaffected.

## 5. Passing of risk

Insofar as no agreement to the contrary has been made the delivery will be ex works in line with the international provisions of INCOTERMS 2010 which form the basis also for all other delivery terms.

## 6. Material defects

We are liable for material defect as follows:

- 6.1 The purchaser agrees with us that the most cost-effective variation will be selected in the case of an additional claim to performance by the purchaser (rectification of defect or additional delivery) insofar as the purchaser does not incur any disadvantages as a result.
- 6.2 Claims in respect of material defects become statute-barred 12 months after the passing of the risk (section 5). This does not apply insofar as longer periods are stipulated according to the Austrian code of Obligations.
- 6.3 The purchaser is required to give notice to us of material defects immediately, though no later than within seven days from the service provision (apparent defects) or the discovery of the defect, in writing. Otherwise the assertion of claims in respect of defects is debarred.
- 6.4 In the case of notices of defect payments by the purchaser may be withheld by the purchaser to an extent appropriately commensurate with the material defects that have occurred. The purchaser may withhold payments only if a notice of defect has been put forward in respect of the justification of which there is no doubt. Should the notice of defect by wrongly issued we shall be entitled require the purchaser to reimburse us for the expenditure incurred.
- 6.5 First of all, we should be provided with an opportunity for additional performance within an appropriate period.
- 6.6 Claims based on defects cannot be lodged in the case of natural wear or damage arising after the passing of the risk as a result of faulty or negligent treatment, excessive stress, and unsuitable means of operation or as a result of particular external influences not anticipated according to the contract. Should improper modifications or maintenance work be carried out by the purchaser or third parties, no claims based on defect can be lodged also in respect of these and the resulting consequences.
- 6.7 Claims by the purchaser in respect of the expenditure for purposes of additional performance, in particular transport, toll, labour and material costs are debarred insofar as the expenditure rises due to the fact that the object of the delivery has subsequently been transported to a location other than the purchaser's establishment unless the transportation corresponds to an appropriate use.
- 6.8 Statutory rights of recourse of the purchaser vis-à-vis ourselves shall arise only insofar as the purchaser has not made any agreements with its customer going beyond claims for defect based on Austria law.
- 6.9 Section 8 shall be applicable to claims for damages. Claims for damages as a result of a material defect going further or other than those dealt with in this section or in section 8 are debarred.

## 7. Information Obligation in case of incidences and quality concerns

- 7.1 All DENDIA GMBH customers and particularly distribution partners of DENDIA GMBH are obliged, to report quality concerns but principally incidences with potential impact – direct or indirect – in life or health of humans immediately and to the attention of DENDIA GMBH Austria!

## 8. General conditions for the return of goods

(With the exception of returns as described in section 6)

- 8.1 Before being able to approve a return we require the following information: the full item number, the quantity, designation, batch number of each package/sales unit. The packaging must be undamaged, unopened and not labelled with any additional labels/identification. We accept only complete repackaged units (no individual packs).
- 8.2 Non-acceptable products are the following: all specially produced products (customer products) or all specially packaged products, all products older than 2 years (according to batch number or delivery note).
- 8.3 The costs of control, quality assurance, possible repackaging and administration are payable in full by the purchaser. We shall deduct 20% of the value of the goods in respect of these costs. All costs in respect of returns, i.e. postage (transportation), insurance, duty, taxes etc are payable in full by the purchaser
- 8.4 The purchaser will be responsible for the return.

## 9. Industrial property rights and copyrights, defects of title

Unless otherwise agreed, we shall be obliged to effect the delivery free from industrial property rights and copyrights of third parties (in the following referred to as property rights) only in the country of the place of delivery. Insofar as a third party lodges justified claims against the purchaser as a result of infringement of property rights by deliveries made by us and utilised in line with the contract, we shall be liable vis-à-vis the purchaser within the period laid down in section 6.2 as follows:

- 9.1 We shall obtain, at our discretion and cost, a right of use for the supplies concerned, modify them in such a way that the property right is not infringed or exchange them. Should this not be possible subject to appropriate conditions, the purchaser shall be entitled to the statutory right of rescission or reduction of purchase price. Our obligation to provide compensation is determined according to section 8.
- 9.2 The above-mentioned obligations arise only insofar as the purchaser notifies us immediately in writing regarding claims asserted by third parties, does not acknowledge that an infringement has taken place and we retain the right to effect all aversion actions and settlement negotiations. Should the purchaser cease utilisation of the supply for purposes of reducing losses or for other important reasons, it is obliged to draw the attention of the third party to the fact that terminating the utilisation is not associated with any recognition of an infringement of property rights.
- 9.3 Claims by the purchaser are debarred insofar as it is responsible for the infringement of property rights.
- 9.4 Claims by the purchaser are also debarred insofar as the infringement of property rights is caused by special demands of the purchaser's, by a use not foreseeable by us or by the supply having been modified or used together with products not supplied by us.
- 9.5 In the case of infringements of property right, the provisions of sections 6.4, 6.5 and 6.9 are applicable correspondingly to the purchaser's claims dealt with in section 7.1.
- 9.6 Further claims of the purchaser or claims other than those dealt with in section 7 vis-à-vis us or persons employed by us in respect of a defect of title are debarred.

## 10. Overall liability

- 10.1 We are liable for damages and compensation for fruitless expenditure in the sense of the Austrian code of Obligations (hereinafter referred to as "damages") as a result of defects of the supply or services or as a result of the infringement of other contractual or non-contractual obligations, in particular as a result of tortious acts only in the case of intent or gross negligence. The above limitation of liability does not apply to loss of life, bodily injury or injury to health, on giving of a guaranteed or accepting a procurement risk, the infringement of essential contractual obligations and in the case of liability according to the product liability law.
- 10.2 Damages due to infringement of essential contractual obligations are limited to compensation for damage which we should have anticipated on conclusion of the contract as a result of circumstances recognisable by us as possible consequence (damage typical of the contract) insofar as there is no intent or gross negligence or as a result of loss of life, bodily injury or injury to health or the giving of a guarantee or a procurement risk.
- 10.3 All limitations on liability apply to the same extent to employees.
- 10.4 The above provisions do not entail a modification of the burden of proof to the detriment of the purchaser.

## 11. Reservation of title

- 11.1 The objects of the delivery (reserved property) remain our property up to the fulfilment of all claims to which we are entitled vis-à-vis the purchaser as a result of the business relationship. Insofar as the value of all security interests to which we are entitled vis-à-vis the purchaser exceeds the level of all secured claims by more than 10% we shall release a corresponding portion of the security interest at the purchaser's request
  - 11.2 During the existence of the reservation of title the purchaser shall not be permitted to effect a pledge or an assignment as security and the resale shall be permitted only to resellers in the ordinary course of business and only subject to the condition that the reseller receives payment from its customer or makes the reservation that ownership passes to the customer only when it has fulfilled its payment obligations.
  - 11.3 In the case of pledges, seizures or other dispositions or interventions by third parties, the purchaser is obliged to notify us immediately to allow us to take legal action. Insofar as the third party is unable to reimburse us for the legal and extra-judicial costs of an action the purchaser shall be liable for the loss incurred by us.
  - 11.4 The purchaser shall be obliged to treat the purchased item with care; in particular, the purchaser shall be obliged to insure it adequately to the value as new at its own expense against damage by fire, water and theft. Insofar as maintenance and inspection operations are required, these will need to be carried out by the purchaser in due course at its own expense.
  - 11.5 In the case of infringement by the purchaser of its obligations, in particular payment delay, we shall be entitled to withdraw from the contract and to take back the goods; the purchaser shall be obliged to surrender them. Taking back of the goods or the assertion of the reservation of title does not require the rescission of the supplier; these actions or pledging of the goods subject to reservation of title by us do not represent a rescission of the contract unless this is expressly declared by us.
  - 11.6 Should the purchaser have sold the purchased item in the ordinary course of business, it passes on to us at this time all outstanding accounts to the amount of the final invoice (including VAT) of our claims arising for the purchaser from the resale vis-à-vis its customers or third parties irrespective of whether the purchased item has been sold on as such or after processing.
- The purchaser remains authorised to collect the receivables under the claim even after the assignment of the claim. This does not affect our authority to collect the receivables under the claim ourselves. However, we undertake to refrain from collecting the receivables for as long as the purchaser complies with its payment obligations by using the proceeds received, does not default in payment and, in particular, no petition for insolvency proceedings is filed or a stoppage of payments arises. However, should this be the case, we shall be able to require the purchaser to inform us regarding the assigned claims and the debtors involved, to provide all the information necessary for collection, to hand over the relevant documents and to inform the debtor (third party) of the assignment.
- 11.7 Processing or transformation of the purchased item by the purchaser will at all times be carried out on our behalf. Should the purchased item be processed together with other objects not belonging to us, we acquire the co-ownership of the new item in proportion of the value of the purchased item (final amount of invoice) to that of the other objects processed, at the time of processing. For the resulting object, the same shall apply to the object created by processing as to the object of sale delivered under reservation. The same applies to the item produced by processing as for the purchased item supplied subject to reservation.
  - 11.8 Should the purchased item be mixed inseparably with other objects not belonging to us, we shall acquire co-ownership of the new item in proportion of the value of the purchased item (final amount of invoice, including VAT) to the other admixed objects, at the time of mixing. Should the mixing operation take place in such a way that the purchaser's item should be regarded as the main item, it is deemed to have been agreed that the purchaser transfers to us co-ownership in proportion. The purchaser preserves the sole ownership or co-ownership which has thus arisen, for us.

## 12. Place of performance, venue, law applicable

- 12.1 Our registered office shall be the place of performance and jurisdiction for all rights and obligations arising from our deliveries and services. We are optionally entitled to sue at the buyer's place of business.
- 12.2 Austrian law shall apply exclusively. The application of the UN Convention on Contracts for the International Sale of Goods (CISG) is excluded.
- 12.3 The data provided by the Buyer shall be stored and processed by electronic data processing, insofar as this is permissible under the Data Protection Act.

Issued: Jan 2021/GL


## DENDIA

© Dendia GmbH, Feldkirch, Austria.  
All rights of distribution, including phototechnical reproduction, reprinting (including excerpts) or storage and recovery in data processing systems of all kinds are reserved and require written permission. Product, color changes, typesetting and printing errors reserved.

Edition: January 2024

Art. No.:108010.5950.0001

*Quality is not a characteristic  
it is a claim!*

**DENDIA**  
DENTAL INSTRUMENTS

**DENDIA GMBH**

Reichsstraße 126  
6800 Feldkirch, Austria

T +43 (0)5522 819 44  
sales@dendiadental.com

WWW.DENDIADENTAL.COM

**Your local dealer:**